


INTRAOCULAR LENS

LUCIS®

HYDROPHILIC ACRYLIC

Incision 2.2 mm



MADE IN **FRANCE**
CE marked

Cristalens

DESIGNATION	TECHNICAL SPECIFICATIONS
Lens type	For implantation in the capsular bag
Optic diameter	6.00 mm
Overall diameter	11.00 mm (from -9.0D to +9.5D) 10.50 mm (from +10.0D to +30.0D)
Design	One piece
Angulation	10°
Material	25% hydrophilic acrylic
Dioptric powers	From -9.0D to +30.0D by 0.5D
Estimated A-Constant (SRK-T)	118.0 Ultrasound biometry 118.5 Interference laser biometry
Suggested Anterior Chamber Depth (ACD)	4.96 mm Ultrasound biometry 5.25 mm Interference laser biometry
Refractive index	1.46
Sterilization	Steam
Recommended incision size	2.2 mm

16/11/2015 DOC_1_PHILIC



Intraocular Lenses
Manufacturer

Cristalens



Cristalens Industrie

4 rue Louis de Broglie - 22300 Lannion - France
Tel: + 33 (0) 2 96 48 92 92 - Fax: + 33 (0) 2 96 48 97 87

www.cristalens.fr