

OPHTHALMIC DIAGNOSTIC STRIPS

In Medical Grade Tyvek Pouches

Schirmer Strips


CE
0434



Sterile Ophthalmic strips for Dry Eye test

Packaging:

Each box contains 100 sterile strips.



Indications:

- Used for measurement of tear production
- Gauged for easy reading

Fluorescein Sodium Ophthalmic Strips

CE
0434




Composition (per strip):

Each strip is impregnated with 1.0 mg of Fluorescein Sodium

Packaging:

Each box contains 100 sterile strips.



Indications:

- For Applanation Tonometry
- For evaluating Contact Lens fitting

High Molecular Fluorescein Sodium Ophthalmic Strips




Composition (per strip):

Each strip is impregnated with 0.5 mg of high molecular Fluorescein Sodium

Packaging:

Each box contains 100 sterile strips.



Indications:

- For Applanation Tonometry
- For evaluating Contact Lens fitting

Lissamine Green Ophthalmic Strips




Composition (per strip):

Each strip is impregnated with 1.5 mg of Lissamine Green

Packaging:

Each box contains 100 sterile strips.



Indications:

- To stain preferentially damaged or devitalized cells
- To indicate any dry patches, as well as any mucus-deficient or damaged Corneal Epithelial cells
- It does not stain healthy, proliferating Ocular surface cells and has minimal effect on cell viability

Rose Bengal Ophthalmic Strips




Composition (per strip):

Each strip is impregnated with 1.5 mg of Rose Bengal


Packaging:

Each box contains 100 sterile strips.



Indication:

- It stains dead or degenerated Epithelial cells of Cornea & Conjunctiva.



BIOTECH
VISION CARE

Manufacturing Facility

555 / 556 / 557, Nr. Shubham Tex-O-Pack,
Khatraj, Dist. Gandhinagar, Gujarat, India.

India Business

210, Sarthik - II, Opp. Rajpath Club,
Sarkhej-Gandhinagar Highway,
Ahmedabad-380054, Gujarat, India.
Ph. / Fax: 079-40038374,32420133
E-mail: sales@biotechvisioncare.com

International Business

401, Sarthik - II, Opp. Rajpath Club,
Sarkhej-Gandhinagar Highway,
Ahmedabad-380054, Gujarat, India.
Ph.: +91-79-26870896/97/98 **Fax:** +91-79-40033767
E-mail: intbusiness@biotechvisioncare.com